









Tissue Viability Service Update/ Covid 19 / Core Ward Product Information

These products below are for standard stock to be available at ward level – this is not an exhaustive list and further information is available on TV webpages on StaffNet.

Product	Indications for use	Contraindications	Application & wear time
Autraman  7x10cm 233182 box of 20 10x20cm 052301 box of 20	Non adherent, polyester mesh wound contact layer. Prevents adherence of granulating tissue onto dressings. Apply as a primary wound contact layer for a wide variety of wounds including infected wounds.	None Listed	Apply directly onto the wound. Effective for up to 7 days (Petrolatum free)
Comfeel Transparent  5x7cm 158140 box of 10 10x10cm 165438 box of 10	Superficial low exuding wounds to debride low levels of slough	Known sensitivities to carboxymethylcellulosegelatin, pectin	Dry skin, remove backing layer and place directly on wound bed, allowing a 3cm border around wound bed. Change when gelling of 80% of dressing has taken place. Remain up to 7 days.
Urgoclean  6 x6cm 214761 box of 10 10x10cm 210510 box of 10 2.5x40cm 214754 box of 5	A sterile, non-adherent, slough-trapping, poly-absorbent fibre dressing with the TLC healing matrix to promote wound healing and enable pain-free dressing changes. The slough trapping fibers (poly-absorbent) bind and trap the slough within the dressing, providing safe and effective desloughing.	Known sensitivities to any components of dressings.	Remove backing film from TLC layer and apply direct to wound bed. Document how many products used to ensure all products removed. Remain 3 – 5 days or dependant on wound exudates.
Urgotul Absorb Border  8 x8cm 211739 box of 10 10x10cm 211722 box of 10 Sacral 214426 box of 10	An absorbent adhesive foam dressing. For all types of exuding wounds.	Known sensitivities of any components of the dressing.	Remove back paper apply directly to wound. Renew according to exudate levels.
Zetuvit Plus  10x10cm 184279 box of 10 20x25cm 184309 box of 10 20x40cm 220595 box of 10	A high capacity absorbent wound pad that holds fluid away from the skin. Moderate to heavily exuding wounds.	Can be applied directly over the wound or as a secondary dressing. Do not use on arterial bleeds or heavily bleeding wounds Discontinue use if allergic to any components of the dressing.	Place white face down on wound surface with at least a 2cm border overlap around the wound area. Do not cut dressings. Wear time dependant on exudates levels. Change if there is staining on the bandage or seepage from the dressing.
Hydroclean plus  5.5cm round HYD770W box of 10 7.5x7.5cm HYD771S box of 10 10x10cm HYD772Y box of 10	Debridement Dressing Hydro responsive dressing used for cleansing absorption and debridement	No Contraindications listed.	Apply directly to wound with film dressing to secure. Can be left in place for up to 3 days between dressing changes dependant on exudate levels.
Prontosan Gel X  50g 193580 each 250g 171743 each	Cleanses and moisturizes chronic and acute skin wounds Prevents the formation of biofilm. Facilitates gentle dressing replacements, Is compatible with commonly used wound dressings.	Known sensitivities of any components of the dressing. Do not use near eyes, inner and middle ear or in combination with soaps, oils cleaning solution	Gel should be applied in thick layer, to cavities and pockets should be filled with gel. Renew application according to wound type. Shelf life 8 weeks post opening.
Film dressing 365  6 x7cm 170944 box of 100 10 x12cm 170951 box of 50 15 x20cm 170968 box of 50	Shower proof film dressing, primary and secondary dressing.	Do not stretch the dressing on application as can cause skin tension.	Dry skin and select suitable size for 2cm overlap around wound edges. Can remain for 7 days.



COVID-19 APPROVED GUIDANCE

OFFICIAL SENSITIVE

Note: This guidance has been fast-tracked for approval for use within NHSGGC

Covid-19 Tissue Viability Core Ward Product Information

This guidance is intended to assist healthcare professionals in the choice of disease-specific treatments.

Clinical judgement should be exercised on the applicability of any guidance, influenced by individual patient characteristics. Clinicians should be mindful of the potential for harmful polypharmacy and increased susceptibility to adverse drug reactions in patients with multiple morbidities or frailty.

If, after discussion with the patient or carer, there are good reasons for not following guidance, it is good practice to record these and communicate them to others involved in the care of the patient.

Version Number:	2
Does this version include changes to clinical advice:	N/A
Date Approved:	23 rd April 2020
Approval Group:	Covid-19 Tactical Group

Important Note:

The version of this document on the Clinical Guideline Directory is the only version that is maintained. Any printed copies should therefore be viewed as 'Uncontrolled' and as such, may not necessarily contain the latest updates and amendments.